

Pensky Marten Flash Point Apparatus (Closed Type)

To determine the flash point of the sample whose flash point is between **40°C to 360°C**. It is widely used to determine the flash point of lubricating oil, bitumen fuel oil, asphalts except cutback bitumen having flash point, biodiesel fuel etc. The apparatus consists of a cast iron stove fitted with a polished brass cover. Brass cup also contains a stirrer for homogenous mixing to maintain the sample temperature. Heater with a suitable energy regulator is provided to the apparatus.

Standard: IP 34, ASTM D 93, IS 1448 (P-21), ISO 2719

Raj Make

Package Includes:

Pensky marten flash point apparatus ,
thermometer range of -5°C to 110°C (IP 15) and
range of 90°C to 370°C (IP 16) and manual.

Ordering Informations:

Ref No.515 Pensky Marten Flash Point App.

HSN Code: 90318000

