

Tri Gloss Meter

Standard: IS 101 (Part 4/ Sec.4), ISO 2813, ASTM D523, ASTM D2457, DIN 67530

Features:

- Measurement range 20°, 60° and 85°.
- Single button for all angles measurement.
- Display 5 readings for measurement data and good for comparison.
- Connect to PC for more extended functions with GQC6 Quality Control Software.
- Large screen display for easy operation.
- Ergonomics design with good appearance.
- Good user communication interface.
- Powerful functions to meet many different requirements.
- Large storage to save over 1000 data.
- Built-in lithium-ion rechargeable battery.

Technical Specifications:

Measuring Angle	: 20°/60°/85°		
Measuring Area	: 20°	60°	85°
	10 x10mm	9 x15mm	5 x36mm
Measuring Range	: 20°	60°	85°
	0~1000GU	0~1000GU	0~160GU
Resolution	: 0.1GU		
Repeatability	: 20°	60°	85°
	± 0.1GU	± 0.2GU	± 0.2%GU
Measuring Time	: 1.0s		
Dimension (L x Wx H)	: 160 x 75 x 90 mm		
Weight	: 350g		
Battery	: 3.7V, 3200mAh Li-ion Battery		
Display	: TFT 3.5 inch, capacitive screen		
Interface	: USB/RS-232		
Storage	: 1000 readings		
Software	: GQC6 Quality Control Software		

Package Includes:

Gloss Meter, Power Adapter, USB Cable,
Standard Calibration Plate and Manual.

Ordering Informations:

Ref No : 225/T Tri Gloss meter

HSN Code: 90275090



Website



Contact

Raj Scientific LLP

Head Office:

201/202 Bhola Bhagwan Industrial Estate, I. B. Patel Road,
Goregaon (E), Mumbai - 400 063.

Work Shop:

119, 136 Shanta Industrial Estate, I. B. Patel Road,
Goregaon (E), Mumbai - 400 063.



: +91-22-2686 3333 / 2685 2115



: +91 98199 04256 / 99677 56248



: info@rajscientificllp.com



: www.rajscientificllp.com