

4 Sided Film Applicator

Applicator are designed for the production of uniform films of paints, adhesive and similar products on plane substrates. They combine the accuracy of fixed applicators with the versatility of multiple clearance/gap choices in one unit. These applicator suitable for use of aqueous acid and alkaline products.

Film width : 4"

Ref.: 113/1 Range 30, 60, 90, 120 Microns

Ref.: 113/2 Range 25, 50, 75, 100 Microns

Ref.: 113/3 Range 50, 100, 150, 200 Microns

Ref.: 113/4 Range 5, 10, 15, 20 Microns

Ref.: 113/5 Range 100, 200, 300, 400 Microns

