

ASTM Solvent Distillation Apparatus

The sample evaporated and condensed under controlled conditions and observations are made of the temperature at various percentage are recovered.

As per IS 1448 (P-18) IP 123 ASTM D 86

Electrically operated single test unit supplied with water bath, measuring cylinder, round bottom flask and Thermometer 10°C to 360°C

Included Accessories

Thermometer

Flash as per IP Standards

Heating coil

Ceramic board with hole (asbestos plate)

Flask 125 ml with M. Cylinders 100ml

Silicon Cork

Ref.: 778

