

Abrasion Tester

The Abrasion Tester is used to test a wide spectrum of materials; The instrument has been referenced in numerous standards and specifications (including plastics, coatings, laminates, leather, paper, ceramics, carpeting, safety glazing, etc.) Two genuine abrasive wheels, which are applied at a specific pressure, are lowered onto the specimen surface. Characteristic rub-wear action is produced by contact of the test specimen against the sliding rotation of the two abrading wheels.

Consist of :-

- Auxiliary Weight : 250 gms. Load (each arm)
- Counter Weight : 50 and 125 gms
- Specimen Adjustment Discs 7 no's
- Calibraes wheel set (similar to CS - 10)
- S-42 Sandpaper Strips Taber USA Make
- Vacuum Unit with suction Hose

