

Digital Muffle Furnace

- Useful for Ashing, Igniting and heat treating of small parts in chemical, industrial fields.
- Outer body made up of heavy gauge M.S./G.I duly powder coated. Inner muffle with high temp. insulating material.
- Maximum Temp. 1000°C / 1200°C
Working temp 930°C / 1130°C
- Temperature Controlled by Digital PID Temp. Controller cum Indicator.
- Heating elements are made of kenthal wire.

	Size	TEMP. (1000°C)	TEMP. (1200°C)
Ref.: 705/1	9" X 4" X 4"	1500 w	2000 w
Ref.: 705/2	10" X 5" X 5"	2000 w	2500 w
Ref.: 705/3	12" X 6" X 6"	3300 w	3500 w
Ref.: 705/4	18" X 9" X 9"	6000 w	7000 w

