

Digital Stormer Viscometer

For measuring viscosity of non newtonian liquid such as paints, resins etc. Rotating paddle driven by adjusting weights to determine viscosity.

- Easy to use.
- Rotation at 200 rpm with audio and visual indication
- Rapid and accurate measurement of viscosity in krebs or centipoise
- Supplied with container and different weights
- Tachometer Version
- Digital RPM Display
- Complies with IS : 101 (Part1/Secs) 1989 or ASTM D 562
- Range 800 - 5220 CPS
 54 - 140 KU

Ref.: 248/1 ASTM D 562

Ref.: 248/2 IS : 101

