

Hot Air Oven

- Inner Chamber made of S.S. and Outer made of Mild Steel duly powder coating
- Gap between the walls filled with glass wool or puf insulation.
- Heating elements are made of high grade chrome plated nichrome wire.
- Temp. Range : 50° C to 250° C
- Accuracy: $\pm 1- 2^{\circ}$ C

Accessories :

- Air Circulation Fan for Uniform Temperature
- Digital PID Temp. Indicator cum Controller
- Digital Timer
- Temp. Range : Up to 300° C / 400° C

Ref.: 701/1 14" x 14" x 14"

Ref.: 701/2 18" x 18" x 18"

Ref.: 701/3 18" x 18" x 24"

Ref.: 701/4 24" x 24" x 24"

Ref.: 701/5 24" x 24" x 36"

