

Hot Plate

Round

Made of cast iron top and body is of thick mild steel painted attractively.

Thermostatic Control Top M. S. / Cast Iron Plate

Size	Rating
8" Dia	1.5 Kw
9" Dia	1.6 Kw
12" Dia	2.0 Kw

Rectangular

Body made of thick mild steel Painted attractively.

The thick iron plate is provided on top.

Thermostatic Control Top M. S. / Cast Iron Plate

SIZE	Rating
10" x 12"	1.5 K.W
10" x 16"	1.5 K.W
12" x 18"	2.0 K.W
12" x 18"	2.0 K.W
18" x 24"	3.0 K.W

