

## Humidity Cabinet

- Double walled, Interior made of stainless Steel exterior made of mild steel duly powder coated.
- Two Stainless steel removable trays.
- Digital temp. cum humidity controller.
- Forced air circulation for uniform temp. In case of cooling models by Hermetically sealed compressor provided.
  
- Temperature Range : 5° C above ambient to 60° C  
10° C to 60° C
- Humidity Range : 40% to 95% RH
- Temperature Accuracy :  $\pm 1^{\circ}$  C
- Humidity Accuracy :  $\pm 3\%$  RH
- Direct RH Controller.
- RH Sensor Range 40 to 100% for Direct RH Controller
- Suitable voltage stabilizer



	<b>Capacity ( Ltrs. )</b>	<b>Inner Chamber size ( cm )</b>
--	---------------------------	----------------------------------

Ref.: 444/1	91	45 x 45 x 45
Ref.: 444/2	175	70 x 50 x 50
Ref.: 444/3	288	90 x 55 x 55
Ref.: 444/4	342	90 x 60 x 60
Ref.: 444/5	441	90 x 70 x 70