

Kinematic Viscometer Bath

As per IP 71, IS 1448 & ASTM D 445-72.

Kinematic Viscometer bath is rectangular shapes Inside tank is made of stainless steel sheet. Two sides toughen glass window provided for full visibility with light arrangement.

Four viscometer to accommodate, double glass with Digital Temperature controller cum indicator. Supplied with viscometer Holder and suction bulb.

Viscometer glass tube canon Fenske direct Flow, reverse flow as per IP and ASTM D standard with constant of bulb.

Size 25, 50, 75, 100, 150, 200, 300, 350, 400, 450, 500, 600

Temperature Range : 5° above ambient to 150°C

Accuracy : ± 1°C

Resolution : 0.1°C

Size : 12" x 9" x 12"

Optional : GMP Model also available

Ref.: 274

