

## Oil Bath

Double walled, outer made of mild Steel duly powder coated  
Inner tank made of Stainless Steel. The bath is provided with  
two heaters of 1 K. W. each rating. Thermo Regulator is  
provided for temperature control. A powerful induction motor  
drives the centrifugal pump.

Temperature 5°C above ambient to 200°C.

Microprocessor PID temperature controller cum indicator  
with timer.

Working Size : 12" x 12" x 8" (L x W x H)

Capacity : 17 Ltr.

Ref.: 408

