

## Penetrometer

For determine semi solid to solid materials, such as grease, waxes, cosmetics, Paint. Using Penetration in the sample by a cone of definite weight (50 or 100) gms fitted to rod. By adjusting arm to touch the surface of the sample, automatic timer is also available.

**Complies with** : ASTM D 937, ASTM D 217, IS 1448, IP 50.

- Centering device
- Dial with brass cone and container

**Ref.: 140/1** Unit with standard accessories and geared arrangement.

**Ref.: 140/2** Unit with standard accessories, Automatic digital timer and geared arrangement.

