

Vibro Shaker

The VIBRO SHAKER is a very useful machine for effective mixing of pastes with liquids, dry powders with solvents, settled pigments with immiscible liquids.

Features:

- Easy to Operate.
- Negligible Maintenance.
- Simultaneous Oscillating and Tilting mechanism.
- Digital Timer with alarm and memory.
- Single or double clamp version oscillation $28^\circ \pm 1^\circ$ angular. Motor $\frac{1}{2}$ HP. Power Supply 230 volts A.C., 50 Hz

Ref.: 810

