



# Density / Viscosity

## Weight Per Litre Cup

A Stainless Steel Precision instrument for determining the specific weight (density) of paints and similar products. This instrument consists of a cylindrical container and cover with a hole for exhaust of excess liquid. Three capacities are available : 25, 50, or 100 cc. These can be supplied with or without a tare weight for weighing on a balance.

Complies with BS 3900, ISO 2811, ASTM D 1475.

Ref.: 120/25      cap. 25 ml

Ref.: 120/50      cap. 50 ml

Ref.: 120/100      cap. 100 ml

