

Wet Abrasion Scrub Tester

Wet Abrasion Scrub Tester, is mainly intended for the determination of the washing and scrubbing resistance as well as cleanability of dispersion paints (Specifically “wall paints”)

Technical Data

- Stoke Rate : 38 ± 2 cycle/min
- Stroke Length : 300 ± 5 mm
- Four Digit Present
- Dimension (W X D X H) : (450 x 350 x 230) mm
- Power Supply : 230 vac

Accessories

- One nylon brush and one sponge brush
- 500 ml washing agent
- Injecting Syringe
- PVC foil, black, fine matt finish size (435 x 165) mm
- Glass Panel

