

Wire Bar Coater

Wire wound rods provided an economical ways to assure the uniformly of thickness of very thin films.

Bar and wire manufactured of stainless steel perfect for very thin films, e.g. foil and paper coating test ideal for flexible materials such as paper, cardboard, test charts, foils, leathers, textiles etc.

Maximum attainable wet film thickness engraved on bar end.

Wet Films Thickness : 4, 10, 15, 25, 40, 50, 70, 80,
100 micron

Film width : 12 or 22 cms.

Impression pad : Small 8.25" x 5"

Rubber on Top Large 12" x 8"

Pad with Glass plate : Medium 10" x 8"
on Top.

Ref.: 115

